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Jen-DuaCem SE

(EN)

Jen-DuaCem SE is recommended for cementing fiber splinting systems, metal, fiberglass, carbon fiber and quartz posts, as well as inlays, onlays and veneers. The substance has improved mechanical properties and radiopacity, in most cases it does not require pre-etching and adhesive preparation. Jen-DuaCem SE is supplied in a two-cylinder syringe with nozzles ensuring automatic and quality mixing of the material as well as high speed and accuracy of application.

●● INTENDED USE

Jen-DuaCem SE is a restorative material intended for a wide range of clinical luting procedures, especially where adequate light curing cannot be guaranteed. The material can be used as:

- Non-light-curing bonding cement for luting indirect restorations.
- Chemical-curing luting composite cement for fixing cast cores.
- For cementation of inlays, onlays and veneers.
- For cementing metal, fiberglass, carbon fiber and quartz pins.
- Luting cement for adhesive bridges like Maryland or Rochette.
- For fixing splinting elements.
- For splinting teeth (with splinting tapes and cords) in cases of illness or injury to periodontal tissues.

●● TECHNICAL INFORMATION

- Polymerization mode: double polymerization.
- Working time: at least 90 seconds at a temperature of 23°C.
- Curing time: about 3-5 minutes at 37°C.
- Light curing: 20 seconds for a layer not more than 2 mm thick or 40 seconds for a layer not more than 4 mm thick.
- Please note that the data given are for a temperature of 23°C and a relative humidity of 50%. These values may vary depending on the conditions and period of storage of the material, temperature, air humidity. Increasing the temperature of the material reduces its operating time.

●● DIRECTIONS FOR USE

DOSAGE AND MIXING:

Before using for the first time, it is important to make sure that the material comes out evenly from both parts of the twin syringe. To do this, squeeze a small amount of material from the twin syringe onto the mixing pad to achieve even paste dosage and ensure appropriate mixing with the mixing pad later on.

Use a syringe with clean, disinfected gloves to avoid contamination. Do not reuse the syringe if it is contaminated. Disinfect the Jen-DuaCem SE twin syringe with a medium level disinfectant after each use.

- Turn the cap 90° anticlockwise, align the notch marks on the cap and the syringe and remove the cap.

- Attach the mixing nozzle, also aligning the notches on the nozzle and the syringe.

- Turn the mixing nozzle 90° clockwise and make sure it is correctly attached to the syringe before applying the material. The material should be squeezed out easily.

- Press gently on the piston to mix and apply the preparation. The two pastes are automatically mixed in their respective volumes. Additional manual mixing of the extruded material before use improves mixing quality!

- To remove the mixing nozzle once the material has been applied, turn it 90° anticlockwise, aligning the notches on the nozzle and syringe, and remove the nozzle. Fit the cap by aligning the notches and turn it 90° clockwise. Keep the two-cylinder syringe tightly closed.

- The mixing nozzle must be disposed of after use. Neither the empty syringe nor the used mixer may be reused.

- Clean used metal tools immediately after use with alcohol or acetone.

APPLICATION TECHNIQUE FOR ALL TYPES OF INDIRECT RESTORATIONS (EXCEPT VENEERS):

1. Remove temporary restoration and temporary cement. Thoroughly clean the preparation area, rinse with distilled water and dry by air without oil and water, but do not overdry. Standard inlay preparations should be slightly rounded and slightly beveled. Avoid tapered edges and undercuts. For onlays, an occlusion reduction of 1.5 - 2 mm is required. The buccal and lingual surfaces should be prepared in such a way that the ceramic covers 2 - 3 mm of enamel and is a continuation of the tooth body. The gingival tissue should be pushed back in such a way as to allow the bonding procedure. The pulp chamber should be protected with a calcium hydroxide pad prior to impression taking.

Note: It is important to thoroughly remove the temporary cement and clean the prepared surface as this is the only way to achieve optimal material adhesion. Zinc oxide-eugenol temporary cements should be removed especially carefully, since phenolic substances (which include eugenol) can slow down the polymerization of the material. For complete confidence, it is recommended to use eugenol-free cements. In addition, contact of the material with substances such as hydrogen peroxide, ethylenediaminetetraacetic acid (EDTA) and/or disinfectants should be avoided, since they can interact with the components of the material and affect its polymerization process and/or adhesion strength.

2. Try on the restoration. Jen-DuaCem SE is used for cementation of well-fitted restorations where a thin layer of cement is required. It is not necessary to fix poorly fitted restorations, in which case it is better to redo the work.

3. During bonding procedures, careful control of oral fluids is necessary to prevent contamination of both the prepared area and the cement. A good result can be achieved by using a cofferdam, gum retraction with threads and cotton rolls.

4. Prepare the orthopedic structure or make sure it has been pretreated according to the manufacturer's instructions.

5. Select the desired shade of cement and prepare the Jen-DuaCem SE Dual Syringe for use according to the instructions above.

6. Apply material to the restoration and/or to the prepared tooth surface, depending on the type of orthopedic structure used. Carefully place the restoration onto the prepared site with light pressure so that excess material flows out. Apply moderate pressure and avoid displacement of the restoration.

Note: Ambient light and high indoor temperatures can shorten the material handling time. It is recommended that the lighting intensity of the dental unit be adjusted during the cementation procedure.

7. Remove excess material by waiting for it to gel within 1-3 minutes after application or after 1-2 seconds of exposure to the polymerizer. Avoid excessive material removal, especially at the restoration-cavity interface. A slight excess of material can be left if necessary. Use a scaler or a dental probe or other suitable instrument to hold the structure in place.

Note: when working on large restorations, remove excess material in a timely manner, especially in hard-to-reach areas (gingival margins, proximal areas). In addition, use a barrier gel to prevent the formation of an oxygen inhibited layer. The gel is applied to the fissure areas after removing excess material, but before the material has polymerized.

8. Cure the material using a curing lamp where available or allow the material to cure chemically. For cementation of all-ceramic or composite restorations, cure Jen-DuaCem SE by directing lamp light through the restoration. If the restoration is cemented in this way, more exposure time is required due to the thickness of the restoration. In any case, the material requires about 5 minutes to cure sufficiently, especially in places that are difficult to reach for the light of a photopolymer lamp. Therefore, if further manipulations are required, it is recommended that after removing excess material, it should be polymerized chemically.

9. After the cement has fully polymerised, grind the edges of the restoration or reshape the occlusal surfaces with fine diamond paste. Cemented restorations must be polished definitively in the areas where corrections have been made, using rubber polishing heads and diamond polishing paste.

10. Check occlusion and functional movements and correct if necessary.

APPLICATION TECHNIQUE FOR VENEERS:

1. Prepare the surface of the tooth for veneers.

2. Maintain control of oral fluids and a dry work area by using rubber dam or alternatively cotton swabs and saliva ejector.

3. Before the installation and final fixation of veneers, it is recommended to fit them in using fitting pastes. Evaluate the correctness of the shape and size of the plate, its color and the need for correction of the proximal contact points. Correct and polish, if necessary, with ceramic or composite polishers. It is not necessary to carry out a thorough occlusion check at this stage to prevent damage to the veneers. After fitting, thoroughly rinse the fitting paste with water, dry with a jet of air without oil and water and etch the veneer surface with phosphoric acid for 15-20 seconds. Rinse and dry.

4. Make sure the veneers have been pre-prepared in the laboratory according to the manufacturer's instructions. Do not re-etch them if they have already been etched in the laboratory.

Note: If this has not been done in the laboratory, the inner surface of veneers should be etched within 4-6 minutes with 4% or within 1-2 minutes with 9% -10% hydrofluoric acid or according to the instructions of the veneer manufacturer. Rinse the etched surface thoroughly with water and dry with air without oil and water. A sign of high-quality processing is the homogeneous, matte and clean surface of the veneers. Apply silane to the inner surfaces of the veneers and let sit for 60 seconds, then dry the surface with a jet of oil- and water-free air. Silane-treated restorations should not be rinsed with water.

5. Prepare the tooth surface. To do this, etch the surface with phosphoric acid according to the manufacturer's instructions (e.g., Phospho-Jen AS), rinse thoroughly with water and dry. It is important not to overdry the etched surface.

6. Apply a bonding agent (such as Jen-Unibond or Jen-Unibond SE) according to the manufacturer's instructions.

7. Select the desired shade of cement and prepare the Jen-DuaCem SE Dual Syringe for use according to the instructions above.

8. Apply the mixed cement directly to the veneer surface and gently adjust to the prepared area with light pressure and remove excess cement with an appropriate tool. Particular attention should be paid to the contact areas and gingival margins.

9. Increase the pressure and, while holding it for a few seconds, use a thin light guide to spot polymerize the material in the area of the front surface for a few seconds. At this stage, it is important not to polymerize the material in the contact and gingival margins.

10. Remove excess material and final cure all surface segments of the restoration for 10-30 seconds.

11. Finish and polish the restoration, check the occlusal contacts and correct if necessary.

APPLICATION TECHNIQUE FOR FIXING METAL, FIBERGLASS, CARBON FIBER AND QUARTZ PINS:

1. Perform depulpation and preparation of the space for inserting the pin, using standard techniques. Make sure the area is thoroughly cleaned, rinsed and dried with paper pins.

Note: In addition, contact of the material with substances such as hydrogen peroxide, ethylenediaminetetraacetic acid (EDTA) and/or disinfectants should be avoided, since they can interact with the components of the material and affect its polymerization process and/or adhesion strength.

2. Select the correct pin size and prepare it according to the manufacturer's instructions.

3. Select the desired shade of cement and prepare the Jen-DuaCem SE Dual Syringe for use according to the instructions above.

4. Insert material into the canal and, if desired, apply a thin layer of material to the pin. Quickly insert the pin into the prepared place, allowing free excess material to come out and remove this excess by waiting for the material to gel within 1-3 minutes after application or after 1-2 seconds of exposure to the polymerizer.

Note: due to the increase in the temperature of the material after its introduction into the oral cavity, the rate of its polymerization increases rapidly. Moreover, ambient light and high indoor temperatures can shorten the material handling time. It is recommended that the lighting intensity of the dental unit be adjusted during the cementation procedure.

5. Cure the material using a curing lamp where available or allow the material to cure chemically. In any case, the material requires about 5 minutes to cure sufficiently, especially in places that are difficult to reach for the light of a photopolymer lamp.

6. Cut the pin and restore the stump.

FIXATION OF CROWNS ON IMPLANTS:

1. Remove temporary restoration and temporary cement. Thoroughly clean the preparation area, rinse with distilled water and dry by air without oil and water, but do not overdry.

2. Clean the implant from temporary cement and other materials in a convenient way and use a diamond or any other bur of your choice to roughen the surface of the implant, rinse it thoroughly and dry it with a stream of air.

3. Make sure the restoration fits well.

4. Select the desired shade of cement and prepare the Jen-DuaCem SE Dual Syringe for use according to the instructions above.

5. Apply mixed material directly to the restoration or implant. The material should cover all surfaces so that after careful installation of the structure in place, its excess flows freely.

6. Remove excess material by waiting for it to gel within 1-3 minutes after application or after 1-2 seconds of exposure to the polymerizer. Use a scaler or a dental probe or other suitable instrument to hold the structure in place.

7. Cure the material using a curing lamp where available or allow the material to cure chemically. In any case, the material requires about 5 minutes to cure sufficiently, especially in places that are difficult to reach for the light of a photopolymer lamp.

8. Finish and polish the restoration, check the occlusal contacts and correct if necessary.

●● PRECAUTIONS AND WARNINGS

• The material contains methacrylates, which may cause allergic reactions. Use medical gloves and avoid contact with skin, mucous membranes and eyes. In case of contact, rinse immediately with plenty of water and seek medical advice.

• The automatic mixing nozzles are exclusively disposable and must be disposed of after use.

• The use of materials containing eugenol may inhibit the polymerisation of the composite.

• Light curing is recommended in a thin layer using a halogen or LED lamp of 1000 mW/cm² or more and a blue light emission spectrum of 400-500 nm.

●● CONTRAINDICATIONS

1. Allergic reactions of patients to the substance are known. In the event of an allergic reaction or the appearance of manifestations of patient sensitivity to the components of the material, it is recommended to stop using it and seek medical help.

2. Direct pulp overlap.

3. The impossibility to ensure perfect dryness of the working field.

●● WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty. The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

●● LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

●● STORAGE AND SHELF LIFE

- The material is used at room temperature and can be stored briefly at a temperature not exceeding 25°C.
- Close the cap tightly immediately after use.
- Long-term storage is recommended at 2 - 8°C to ensure maximum shelf life.
- Shelf life: 2 years
- DO NOT FREEZE!

●● RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

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