Fiber dental materials are specially glass ceramic fibers formed in ropes and tapes that can be attached to a tooth tissue with the help of the bonding systems used in your daily dental practice.

## •• INDICATIONS

Flexible enough to bend and almost invisible after coating with resin, the ropes and tapes combine the ability to fine chemical bonding with ease of use, and provide superior aesthetics and durability at implementation of the following clinical procedures:

- Periodontal splinting of anterior and lateral teeth.
- Immediate replacement of lost or missing teeth.
- Stabilization and immobilization of reimplanted or damaged teeth.
- Temporary reinforcement of implant-based bridges.
- Postorthodontic retention.
- Manufacture of temporary and long-term adhesive bridges.

## •• DIRECTIONS FOR USE

To preserve the weaving integrity the ends of the fiberglass material are impregnated with a resin and polymerized by the manufacturer. To properly cut the material prior to use, follow these tips:

1. Place a small amount of adhesive (area with length of 3-4 mm) directly after a portion polymerized at the end.
2. Cure with light.
3. Cut the material by any sharp scissors straight on the hardened area. Now the end of the material will not rumple.
4. Measure the length of material required for this clinical procedure.
5. Place a small amount of adhesive (3-4 mm) on the material in the place specified during the measuring.
6. Carry out polymerization.
7. Cut the material by any sharp scissors straight on the hardened area.
8. Now you have a piece of the desired length with the ends that will not rumple into separate fibers!

## •• BASIC CONCEPTS

- There are used flexible ropes and tapes that are easily bonded to virtually all dental composite resins. The surface of the tooth enamel or dentin is treated with the bonding system, which is now more commonly used in your clinic. No need to change an adhesive system used by you.
- The Jen-Fiber material should be carefully removed from the packaging. Quickly close the cover and avoid contamination of the remaining ropes or tapes. Note that the ends of Jen-Fiber materials are protected. This prevents fiberization of the special weave for intraoral use.
- To prevent fiberization, do not cut the Jen-Fiber material until it is coated with polymerized resin!
- After the material is coated by resin, cut off and discard approximately 3 mm from the tip end to prevent penetration of the material portion polymerized by manufacturer to the final structure.

- Durable flexible element can now be installed in any desired configuration and stowed in a groove etched in a tooth treated with an adhesive composition. Then, to complete the coating of the Jen-Fiber material, use flowable composite resin that may be used by you at your clinic to work on the front and back teeth. We recommend using a flowable composite materials, for example, Jen-LC Flow.
- When you need to replace a tooth, you can construct the bridge body using the fiber-filled reinforced composite, which is to be covered with a layer of microhybrid or nanocomposite for improved aesthetics.
- For the final lining of the bridge body we recommend using of composite materials such as Jen-Radiance or Jen-Favorite LC.

## •• SAFETY WORK

When working with the material, use appropriate protective eyewear, clothing, mask and gloves. Safety glasses are recommended for patients.

## •• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty.

The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

## •• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

## •• STORAGE AND SHELF LIFE

- No special restrictions.
- The shelf life is not affected.

## •• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

## •• PACKING

### Jen-Fiber Tape (2 mm):

- 3 tapes of 9 cm each (total length 27 cm), width 2 mm.

### Jen-Fiber Rope (1,5 mm):

- 3 cords of 9 cm each (total length 27 cm), diameter 1,5 mm.

### Jen-Fiber Rope (3 mm):

- 3 cords of 9 cm each (total length 27 cm), diameter 3 mm.

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EN Jen-Fiber Rope - fiberglass flexible rope for reinforcing the restorations

Jen-Fiber Tape - fiberglass tape for teeth splinting and restorations reinforcing

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