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Jen-Radiance EN

Universal submicron hybrid composite that achieves high aesthetic results using only a small range of shades.

•• INDICATIONS

- Highly esthetic restorations of I, II, III, IV and V cavities;
- Onlay and Inlay;
- “Sandwich” technique;
- Fiber Splinting;
- Crown Build Up.

•• INSTRUCTION

Cavities are prepared in the conventional manner. In very deep or large restorations and in cases of near pulp exposure the use of a calcium hydroxide liner is recommended e.g. Jen-Line LCF.

•• RESTORATING OF ANTERIOR GROUPS OF TEETH

The use of a mylar strip is recommended for interproximal separation and as a matrix. Etch the enamel surrounding the cavity including undercuts with enamel conditioner [Phospho-Jen AS Gel] for 15 – 20 sec., wash and dry. Apply Bonding Agent [Jen-Unibond is recommended] over the enamel margins and cure for 10 seconds. Place restorative material with a tapping movement and cure. In very deep cavities more than one layer may be required. Shape with appropriate finishing diamond or multi-fluted carbide burs, contour proximal surfaces with finishing strips. Check occlusion with a thin articulating paper, carefully adjust occlusion by removing material with a fine polishing diamond or stone. Finish with discs and rubber cups.

•• RESTORATING OF POSTERIOR GROUPS OF TEETH

Place a mylar strip or matrix band for interproximal separation and fix it with wedges. Adapt the band to seal the gingival area to avoid overhangs. Note: the matrix may be placed following the enamel etching if preferred. Apply Bonding Agent [Jen-Unibond is recommended] over the enamel margins and cure for 10 seconds. To aid in adaptation, the first 1,5 mm layer of Jen-Radiance may be placed and adapted to the proximal box. A condensing instrument can be used to adapt the material to all of the internal cavity aspects and cure it. The next portions of Jen-Radiance may be placed in layers less than 2,5 mm. Cure each increment separately. Shape with appropriate finishing diamond or multi-fluted carbide burs, contour proximal surfaces with finishing strips. **Attention!** Jen-Radiance Opaque shades should be applied in layers no more than 1 mm thick.

•• CAUTIONS

Jen-Radiance contains methacrylic resins that may cause an allergic reaction by skin contact in certain individuals. Avoid long or repeating contact of not-polymerized material with the skin (allergic contact dermatitis is possible), soft tissue of the oral cavity and eyes.

In case of contact, immediately wash thoroughly a place of contact with water and soap. If there is rash or other signs of allergic reactions on the skin stop use of the material and ask for medical care.

Use of this product in patients with known acrylic allergies is to be avoided.

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